

A6 22/6/2017 1 of 18

Face and body Analysis Natural Computer Interaction

REST API iMX6 - PC (Controller)







A6 22/6/2017 2 of 18

Contents

1	Overview	2
2	General comments	3
3	Return Values for all REST API calls	3
4	Use Cases / Context	3
	4.1 Initialize	
	4.2 Enrollment	
	4.3 Profile handling	4
	4.4 Use case handling	4
5	iMX6 Server	4
	5.1 Initialize Server	4
	5.2 Set log level	4
	5.3 Enum user profiles	5
	5.4 Create user profile	5
	5.5 Activate user profile	6
	5.6 Delete user profile	6
	5.7 Enum enrollment modules	6
	5.8 Start enrollment of a module	7
	5.9 Enum use case modules	
	5.10 Start use case	10
	5.11 Stop use case	11
	5.12 Setup gaze world model	
	5.13 Define gaze regions	12
	5.14Define gaze + gesture events	13
6	PC Server (server running on the client for callbacks)	14
	6.1 Logging / Warning / Errors	14
	6.2 Use Case Event	
	6.3 Enrollment Data Event	

1 Overview

This document describes the REST API between the iMX6 and the PC.

Section 4 describes which interfaces are used in which context / use case, while 5 and 6 describes the technical interface in detail.

In this document, the iMX6 Hardware will have the IP address 192.168.0.10:8000 while the PC hardware will have the IP address 192.168.0.20:9000.



Rev A6 Date 22/6/2017 Page 3 of 18

2 General comments

- the *Action* tag in each of the transfered JSON structures is technically not needed, but helps when debugging
- every time an interface is called a TRACE event (6.1) should be send back to for debugging purposes.

3 Return Values for all REST API calls

The REST API supports these return value:

200	$\cap K$	request successful
200	UK	request successiui

403 Forbidden page is called, but the server was not initialized before

400 Bad request request has invalid input parameters

404 Not found called page not found

501 Not Implemented called page found, but function has to be implemented

503 Service Unavailable call can not be completed successfully because of internal issues

- more details optionally in returned JSON structure

Each API calls may return these values:

Result Code	Result	Additional Comments
200	call successful	
403	server was not initialized	
400	call failed	
	- further information via LoggingEvent	
501	during development phase	

Additional values are listed separately at each function.

4 Use Cases / Context

This section describes very roughly, which commands are used in with use case / context.

4.1 Initialize

$PC \Rightarrow iMX6$	Initialize	(section 5.1 page 4)
$PC \Rightarrow iMX6$	SetLogLevel	(section 5.2 page 4)
$iMX6 \Rightarrow PC$	LoggingEvent	(section 6.1 page 14)
$PC \Rightarrow iMX6$	SetupGazeWorldModel	(section 5.12 page 11)
$PC \Rightarrow iMX6$	DefineGazeRegions	(section 5.13 page 12)
$PC \Rightarrow iMX6$	DefineGazeGestureEvents	(section 5.14 page 13)

4.2 Enrollment

$PC \Rightarrow iMX6$	EnumEnrollmentModules	(section 5.7 page 6)
$PC \Rightarrow iMX6$	StartEnrollment	(section 5.8 page 7)
$PC \Rightarrow iMX6$	EnrollmentDataEvent	(section 6.3 page 18)





Rev A6 Date 22/6/2017 Page 4 of 18

4.3 Profile handling

$PC \Rightarrow iMX6$	EnumEnrollmentModules	(section 5.3 page 5)
$PC \Rightarrow iMX6$	CreateUserProfile	(section 5.4 page 5)
$PC \Rightarrow iMX6$	ActivateUserProfile	(section 5.5 page 6)
$PC \Rightarrow iMX6$	DeleteUserProfile	(section 5.6 page 6)

4.4 Use case handling

$PC \Rightarrow iMX6$	EnumUseCaseModules	(section 5.9 page 9)
$PC \Rightarrow iMX6$	StartUseCase	(section 5.10 page 10)
$PC \Rightarrow iMX6$	StopUseCase	(section 5.11 page 11)
$PC \Rightarrow iMX6$	UseCaseEvent	(section 6.2 page 14)

5 iMX6 Server

5.1 Initialize Server

5.1.1 IN

URL: http://192.168.0.10:8000/Initialize/Body:

```
1 {"Action": "InitializeServer",
2  "CallbackServer": "192.168.0.20:9000"
3  "OnLoggingEventURL": "/LoggingEvent/"
4  "OnEnrollmentEventURL": "/EnrollmentDataEvent/"
5 }
```

OnLoggingEventURL is the call back URL for logging events on the server running on the client side (see section 6.1).

OnEnrollmentEventURL is the call back URL for the steps during the enrollment phase (see section 6.3).

5.1.2 OUT

no out message - only result code

5.2 Set log level

5.2.1 IN

URL: http://192.168.0.10:8000/SetMinLogLevel/Body:

```
{"Action" : "SetLogLevel",
   "MinLogLevel" : "<Minimal log level>"
}
```





A6 22/6/2017 5 of 18

Valid values for *MinLogLevel* are the TRACE, DEBUG, INFO, WARN, ERROR, FATAL (like in NLOG)

5.2.2 OUT

no out message - only result code

5.3 Enum user profiles

5.3.1 IN

URL: http://192.168.0.10:8000/EnumUserProfiles/
Body: currently optional - if more data has to be passed, it must be added here

```
{"Action": "EnumUserProfiles",
}
```

5.3.2 OUT

```
{"Action": "EnumUserProfilesResult",

"ProfileCount": 3,

"Profiles": [

"<Profile ID1>",

"<Profile ID2>",

"<Profile ID3>"

]

8 }
```

5.4 Create user profile

5.4.1 IN

URL: http://192.168.0.10:8000/CreateUserProfile/Body:

```
{"Action" : "CreateUserProfile",
    "ClientAppUserID" : "<client user ID>"
}
```

<cli>tuser ID> is an arbitrary string that is unique on client side

5.4.2 OUT

```
{"Action": "CreateUserProfileResult",
    "ClientAppUserID" : "<client app user ID>",
    "FANCIUserID" : "<FANCI user ID>"
```





A6 22/6/2017 6 of 18

4

<FANCI user ID> is the ID that is used to store the data of the new user on the IMX6.

5.5 Activate user profile

5.5.1 IN

URL: http://192.168.0.10:8000/ActivateUserProfile/Body:

```
{"Action" : "ActivateUserProfile",
    "FANCIUserID" : "<FANCI user ID>"
}
```

5.5.2 OUT

no out message - only result code

Result Code	Result	Additional Comments
200	Profile is activated	
400	Profile not activated - FANCIUserID unknown	

5.6 Delete user profile

5.6.1 IN

URL: http://192.168.0.10:8000/DeleteUserProfile/Body

```
1 {"Action" : "DeleteUserProfile",
2  "FANCIUserID" : "<FANCI user ID>"
3 }
```

5.6.2 OUT

no out message - only result code

Result Code	Result	Additional Comments
200	Profile is activated	
400	Profile not deleted - FANCIUserID unknown	

5.7 Enum enrollment modules

5.7.1 IN

URL: http://192.168.0.10:8000/EnumEnrollmentModules/

Body: - currently optional - if more data has to be passed, it must be added here

```
1 {"Action": "EnumEnrollmentModules",
2 }
```





A6 22/6/2017 7 of 18

5.7.2 OUT

```
{"Action": "EnumEnrollmentModulesResult",
1
    "ModulesCount" : 3,
2
    "Modules" : [
3
4
        "Module" : "Gaze",
5
        "StartEnrollmentUrl" : "http://192.168.0.10:8000/enroll_gaze"
6
        "Status" : "OK",
7
8
      },
9
        "Module" : "Voice",
10
        "StartEnrollmentUrl": "http://192.168.0.10:8000/enroll_voice"
11
        "Status" : "ERROR",
12
13
      },
14
        "Module" : "Face",
15
        "StartEnrollmentUrl": "http://192.168.0.10:8000/enroll_face"
16
        "Status" : "0K",
17
18
19
    1
20
```

Supported options for *Module* depend on the iMX6 hardware

Supported options for *Status* are: OK, ERROR. Enrollment of modules in Status ERROR can not be started. The Face module is doing enrollment for both Authentication and Emotion.

5.8 Start enrollment of a module

5.8.1 IN

URL: http://192.168.0.10:8000/StartEnrollmentUrl/ - depends on the enrollment module - from EnumEnrollmentModules

Body: Content - currently optional - if more data has to be passed, it must be added here

```
1 {"Action": "StartEnrollment",
2    "Module": "<module name>",
3    "Parameter1": "<value1>",
4    "Parameter2": "<value2",...
5 }</pre>
```

Valid values for *<module name>* are the modules returned in *EnumEnrollmentModules*. *ParmeterX* and the corresponding value *<valueX>* are module specific.

5.8.2 OUT

```
1 {"Action": "StartEnrollmentResult",
```





Rev A6 Date 22/6/2017 Page 8 of 18

```
"Module": "<module name>",
"Status": "<status>",
"Message": "<module specific message>",
"Timestamp": "<timestamp>
6 }
```

Supported options for <module name> are the modules returned in EnumEnrollmentModules. Supported options for <status> are: OK, ERROR.

Supported options for <module specific message> are specific to the module.

<ti>stimestamp> is the timestamp in milliseconds of the iMx6 server. If the status is ERROR, then the call has returned 503.

If the status is OK, then enrollment has started, and the call has returned 200.

Face For Face enrollment, the specific input parameters are the following:

- FANCIUserID the id of the user to enroll.
- *DurationSeconds* the amount of time used for enrollment. If not specified it is default to 5 seconds.

The output *<module specific message>* is one of the following:

- SUCCESS if the Status is OK.
- SERVER NOT INITIALIZED if the Status is ERROR because the server was not initialized.
- SERVER ERROR if the Status is ERROR because the server had an error.
- MISSING_RequiredParameter if the Status is ERROR because the RequiredParameter was not set.
- WRONG_RequiredParameter if the Status is ERROR because the RequiredParameter was not valid.
- NO FACE DETECTED if the Status is ERROR because no face was detected.
- NOT_ENOUGH_FACE_DETECTED if the Status is ERROR because not enough faces were detected.
- MULTIPLE_FACES_DETECTED if the Status is ERROR because multiple faces were detected at the same time.
- *INTERNAL_ERROR* if the *Status* is *ERROR* because one of the authentication or emotion core library failed.
- *CAMERA_BUSY* if the *Status* is *ERROR* because the camera is already used by some other process.
- CAMERA FAILED if the Status is ERROR because the camera could not be opened.
- UNKNOWN ERROR if the Status is ERROR because there was an unknown error.





A6 22/6/2017 9 of 18

5.9 Enum use case modules

5.9.1 IN

URL: http://192.168.0.10:8000/EnumUseCaseModule/

Body: Content - currently optional - if more data has to be passed, it must be added here

```
1 {"Action": "EnumUseCaseModules",
2 }
```

5.9.2 OUT

```
{"Action": "EnumUseCaseModulesResult",
1
2
     "UseCaseCount" : 4,
     "UseCases" : [
3
4
          "UseCase" : "GestureDetection",
5
          "Status" : "OK",
6
          "Enabled" : "true",
7
          "StartUrl" : "http://192.168.0.10:8000/start_gesturedetection"
8
          "StopUrl" : "http://192.168.0.10:8000/stop_gesturedetection"
9
          "ModuleDetails":
10
11
               "SupportedGestures": ["WipeLeft", "WipeRight", "WipeUp"],
12
               "SupportedPostures" : ["ThumpUp", "ThumpDown"],
13
               "..." <mark>:</mark> "..."
14
           }
15
       },
16
17
          "Module" : "GazeTracking",
18
          "Status" : "OK",
19
          "Enabled" : "true",
20
          "StartUrl" : "http://192.168.0.10:8000/start_gazetracking"
21
          "StopUrl" : "http://192.168.0.10:8000/stop_gazetracking"
22
          "ModuleDetails":
23
24
              "TrackerURLs" : [
25
                  "tobii-ttp://is4be-030105141715",
26
                  "tobii-ttp://is4be-030105042734"
27
28
                  ],
              "..." : "..."
29
30
       },
31
32
          "Module" : "Emotion",
33
          "Status" : "OK",
34
```

Rev A6 Date 22/6/2017 Page 10 of 18

```
"Enabled" : "true",
35
          "ModuleDetails":
36
37
              "..." <mark>:</mark> "..."
38
39
       },
40
41
          "Module" : "Authentication",
42
          "Status" : "OK",
43
          "Enabled" : "false",
44
          "StartUrl": "http://192.168.0.10:8000/start_authentication"
45
          "StopUrl": "http://192.168.0.10:8000/stop_authentication"
46
          "ModuleDetails" :
47
48
              "..." : "..."
49
50
       },
51
52
          "Module": "EmotionDetection",
53
          "Status": "ERROR",
54
          "Enabled" : "false",
55
          "StartUrl" : "http://192.168.0.10:8000/start_emotiondetection"
56
          "StopUrl" : "http://192.168.0.10:8000/stop_emotiondetection"
57
          "ModuleDetails":
58
59
               ..." : "..."
60
            }
61
       }
62
63
64
   }
```

Supported options for Module depend on the iMX6 hardware

Supported options for Status are: OK, ERROR.

Supported options for Enabled are: false, true.

The structure of *ModuleDetails* is sensor dependent - this is not JSON conform, but we will deal with that.

5.10 Start use case

5.10.1 IN

URL: http://192.168.0.10:8000/StartUseCaseUrl/ - depends on the use case - from Enum use case modules

Body: Content





Rev A6 Date 22/6/2017 Page 11 of 18

5.10.2 OUT

no out message - only result code

5.11 Stop use case

5.11.1 IN

URL: http://192.168.0.10:8000/StopUseCaseUrl/ - depends on the use case - from Enum use case modules

Body: Content

```
1 {"Action": "StopUseCase",
2 }
```

5.11.2 OUT

no out message - only result code

5.12 Setup gaze world model

5.12.1 IN

URL: http://192.168.0.10:8000/SetupGazeWorldModel/

POST parameter: Content

```
{"Action" : "SetupGazeWorldModel",
1
    "TrackerCount" : 2
2
    "Trackers" : [
3
4
      {
        "active" : "true",
5
        "tracker_url": "tobii-ttp://is4be-030105141715",
6
        "tracker_to_monitor_x_offset_mm" : 0,
7
        "monitor_width_mm" : 292.6,
8
        "monitor_height_mm" : 109.7,
9
        "monitor_width_px" : 1024,
10
        "monitor_height_px" : 768,
11
        "tracker_angle_y_degree" : 0,
12
        "position" : "center"
13
      },
14
15
        "active" : "true",
16
        "tracker_url": "tobii-ttp://is4be-030105042734",
17
        "tracker_to_monitor_x_offset_mm" : -132.5,
18
```





Rev A6 Date 22/6/2017 Page 12 of 18

```
"monitor_width_mm" : 264.5,
"monitor_height_mm" : 165.2,
"monitor_width_px" : 800,
"monitor_height_px" : 600,
"tracker_angle_y_degree" : -10,
"position" : "right"
}
```

5.12.2 OUT

no out message - only result code

5.13 Define gaze regions

5.13.1 IN

URL: http://192.168.0.10:8000/DefineGazeRegions/Body:

```
{"Action" : "DefineGazeRegions",
1
    "InteractorCount" : 2,
    "Interactors" : [
3
4
        "InteractorID" : <InteractorID>,
5
        "ScreenPosition" : <ScreenPosition>,
6
7
        "RectShapeRegionCount" : 2,
        "RectShapeRegions" : [
8
9
             "X" : 10, "Y" : 10, "W" : 10, "H" : 10,
10
             "Weight" : 64
11
          },
12
          {
13
             "X" : 30, "Y" : 30, "W" : 10, "H" : 10,
14
             "Weight" : 64,
15
          }
16
17
         ],
        "CircleShapeRegionCount" : 2,
18
        "CircleShapeRegions" : [
19
20
                  "X": 40, "Y": 40, "R": 10,
21
                  "Weight": 64
22
23
                },
24
                  "X" : 50, "Y" : 50, "R" : 30,
25
```





Rev A6 Date 22/6/2017 Page 13 of 18

```
"Weight" : 64,
26
                }
27
28
              ],
         "Z0rder" : 100
29
30
31
         "InteractorID" : <InteractorID>,
32
         "RectShapeRegionsCount" :1,
33
         "RectShapeRegions" : [
34
           {
35
             "X" : 10, "Y" : 10, "W" : 10, "H" : 10,
36
             "Weight" : 64
37
           }
38
         ],
39
         "CircleShapeRegionsCount": 0,
40
         "CircleShapeRegions" : [],
41
         "Z0rder" : 50
42
44
45
```

A call of DefineGazeRegions always defined updates the complete gaze region set.

The *InteractorID* is defined by the client application.

The *ScreenPosition* is defined in the gaze world model.

The *Weight* is a score that influences how likely a region is to be selected. This is useful if two regions are close to each other, to avoid "flickering" between the two.

5.13.2 OUT

no out message - only result code

5.14 Define gaze + gesture events

5.14.1 IN

URL: http://192.168.0.10:8000/DefineGazeGestureEvents/Body





Rev A6 Date 22/6/2017 Page 14 of 18

```
9
         "EventID" : "<EventID2>",
10
        "InteractorID" : "Button2",
         "Gesture" : "HighFive",
12
      },
13
14
         "EventID" : "<EventID3>"
15
        "InteractorID" : "Button1",
16
         "Gesture" : "WipeLeft",
17
      },
18
19
20
```

5.14.2 OUT

no out message - only result code

6 PC Server (server running on the client for callbacks)

6.1 Logging / Warning / Errors

6.1.1 IN

URL: http://192.168.0.20:9000/LoggingEvent/Body:

```
1 {"Action": "LoggingEvent",
2 "Level": "<log level>",
3 "Message": "<message>",
4 "Description": "<description>"
5 }
```

Valid values for $< log\ level>$ are the TRACE, DEBUG, INFO, WARN, ERROR, FATAL (like in NLOG) < message> is a short description of the logging event.

<description> is an optional detailed description of the logging event.

6.1.2 OUT

no out message - only result code

6.2 Use Case Event

6.2.1 IN

URL: http://192.168.0.20:9000/UseCaseEventUrl/

The URL is the CallbackUrl given when starting the use case (section 5.10). Body:





Rev A6 Date 22/6/2017 Page 15 of 18

```
{"Action": "UseCaseEvent",

"Module": "<module name>",

"Timestamp": "<timestamp>",

"EventData" :

{
    <event data structure - depends on <module name>>

}

}
```

Valid values for < module name > are the Module values returned by EnumUseCaseModules (section 5.9)

The value of *EventData* depends on the sensor module which triggered the event. *timestamp* is the timestamp of the iMX6 server.

Authentication

For the authentication use case, five events can be generated. The json returned is one of the following:

```
1 {"Action": "UseCaseEvent",
2  "Module": "Authentication",
3  "Timestamp": "<timestamp>",
4  "EventData" : {
5     "Name" : "Engagement"
6   }
7 }
```

or

```
1 {"Action": "UseCaseEvent",
2  "Module": "Authentication",
3  "Timestamp": "<timestamp>",
4  "EventData" : {
5     "Name" : "NewUserID",
6     "User" : "<FANCIUserID>",
7     "Score" : "<AuthenticationScore>"
8  }
9 }
```

or

```
1 {"Action": "UseCaseEvent",
2  "Module": "Authentication",
3  "Timestamp": "<timestamp>",
4  "EventData" : {
5     "Name" : "NoFace"
6     }
7 }
```





Rev A6 Date 22/6/2017 Page 16 of 18

or

```
1 {"Action": "UseCaseEvent",
2  "Module": "Authentication",
3  "Timestamp": "<timestamp>",
4  "EventData" : {
5     "Name" : "Impostor",
6     "Score" : "<AuthenticationScore>"
7   }
8 }
```

or

```
{"Action": "UseCaseEvent",
1
   "Module": "Authentication",
2
   "Timestamp": "<timestamp>",
3
   "EventData" : {
4
     "Name" : "SpoofingAttempt",
5
     "User" : "<FANCIUserID>",
6
     "Score": "<AuthenticationScore>"
7
8
  }
```

where FANCIUserID is the ID of the user who was authenticated, or whom someone tried to spoof; and AuthenticationScore is a float value in [-1.0, 1.0] is the score output.

Mood

```
1 {"Action": "UseCaseEvent",
2  "Module": "MoodChange",
3  "Timestamp": "<timestamp>",
4  "EventData" : {
5     "Name" : "<NewMoodName>"
6   }
7 }
```

where NewMoodName is one of NewMoodAnger, NewMoodDisgust, NewMoodFear, NewMoodFear, NewMoodSadness or NewMoodSurprise.

Gesture

```
1 {"Action": "UseCaseEvent",
2  "Module": "Gesture",
3  "Timestamp": "<timestamp>",
4  "EventData" : {
5     "Name" : "NewHandGesture",
```





Rev A6
Date 22/6/2017
Page 17 of 18

```
"Gesture" : "<GestureOrPostureName>",
     }

8 }
```

where *GestureOrPostureName* is one of *BIG_5*, *THUMB_UP*, *V_LETTER*, *FIRST_3*, *FINGERS4*, *PHONE*, *STOP*, *WAVE*, *SWIPE_DOWN*, *SWIPE_UP*, *SWIPE_LEFT* or *SWIPE_RIGHT*. These postures are described in SoftKinetic Close Interaction Library.


Gaze

```
{"Action": "UseCaseEvent",
1
    "Module": "Gaze = ",
2
    "Timestamp": "<timestamp>",
3
    "EventData" : {
4
       "Name" : "GazeData",
5
       "Interactor" : "<InteractorID>",
6
       "Event" : "Enter",
7
       "InteractorStackCount" : 3,
8
       "InteractorStack" : [
9
         "<InteractorID>",
10
         "<InteractorIDBelow1>",
11
         "<InteractorIDBelow2>"
12
13
14
15
```

Valid values for <InteractorID> and InteractorStack are the values uses in DefineGazeRegions (section 5.13)

Valid values for *Event* are the *Enter* and *Leave*

InteractorStack will return the stack of all interactors in the order of the zOrder.



p.e.

- if you gaze at A the value of Interactor should be *Popup* and the value of InteractorStack should be *Popup*, *TabLeft*, *Screen1*
- if you gaze at A the value of Interactor should be *Popup* and the value of InteractorStack should be *Popup*, *TabRight*, *Screen1*





Rev A6 Date 22/6/2017 Page 18 of 18

Gaze-Gesture-Fusion

```
{"Action": "UseCaseEvent",

"Module": "GestureGesture",

"Timestamp": "<timestamp>",

"EventData" : {

"Name" : "GazeGestureData",

"EventID" : "<EventID>",

}
```

6.2.2 OUT

no out message - only result code

6.3 Enrollment Data Event

6.3.1 IN

URL: http://192.168.0.20:9000/EnrollmentDataEvent/

GET parameter: *Module* = < module name > Body:

```
1 {"Action": "EnrollmentDataEvent",
2    "Module": "<module name>",
3    "Progress" : <Progress as integer value 0..100>,
4    "Step" : "<Tag of the enrollment step>"
```

The valid values for *Step* depend on the value of *Module*:

Gaze

5 }

- LookAtGazeCalibrationPostition1
- LookAtGazeCalibrationPostition2
- **–**
- LookAtGazeCalibrationPostitionX

The gaze calibration point coordinates, numbering and amount are entered in the gaze world model.

6.3.2 OUT

no out message - only result code